

# Paper A — data handoff SOP (one page)

As at: 10 Aug 2026 · Pack 1.1.63  
For: Pedro / Ed / John · Beacon feasibility n≈5  
Protocol: PAPER\_A\_FEASIBILITY\_PROTOCOL · kit ORALABLE\_RESEARCH\_KIT · Day-1 ED\_PEDRO\_QUICK\_START

**One-liner:** After each arm, send John the files below (AirDrop / Files / Mail). Do **not** rely on Share → continuous CSV alone for multi-hour arms — that export is only ~3 min of RAM.

## Naming

Oralable\_<ARM>\_<Who>\_<YYYYMMDD>\_<kitID>

Examples: Oralable\_CORE\_Pedro\_20260820\_RK01 · Oralable\_ARMP\_Pedro\_20260821\_RK01 ·  
Oralable\_DUALA\_Mac\_20260822\_RK02

Send: files + one CRF row (feasibility §4). Note FW string (1.0.70) and app (4.3.3).

## Which button / which file (by arm)

| Arm                   | Duration        | In the Oralable app — do this                                                                                                                  | Send John                                                                   | Do not use as sole handoff                              |
|-----------------------|-----------------|------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------|---------------------------------------------------------|
| Core                  | ≥5–10 min       | Connect → temple → HR/SpO <sub>2</sub> stable → <b>Share</b> → <b>Save CSV to Files</b> (or Share sheet). Optional: screenshot vitals card     | CSV + CRF                                                                   | —                                                       |
| Arm P                 | 1–2 h           | Wear continuous. Keep phone nearby. End → <b>Share</b> → <b>Export PDF — Clinical Temporalis Report</b> . Also save any session CSV if offered | <b>PDF</b> (primary) + CRF + MAD fields + AcuPebble notes if any            | Continuous CSV alone (~3 min RAM)                       |
| Arm E/J               | ≥6 h (goal 8 h) | Overnight / long wear. Morning → <b>Clinical Temporalis PDF</b> + hypnogram on Dashboard/Share                                                 | <b>PDF</b> + CRF. Prefer wear ≥6 h before calling it overnight <i>N</i>     | Continuous CSV alone; 1–2 h Arm P is <b>not</b> Arm E/J |
| Dual A (preferred)    | ~6 min          | <b>Mac:</b> run_dual_protocol_a_session.py → then align_anr_oralable_concordance.py                                                            | Pair logs + plots/concordance/<session>/ (NEST.md, overlay) + CRF Dual A ID | Claiming nest as AHI/ODI                                |
| Dual A (iOS optional) | ~6 min          | Settings → tap About 7x → Developer Settings → <b>Dual Protocol A ON</b> → Dashboard → Dual Protocol A → run → <b>Share export pack</b>        | TEMPORALIS_RAW_* + ANR_EMG_* + DUAL_PAIR_* (John runs Mac align)            | Leaving Dual A ON for App Store / Day-1 vitals          |

Sleep is the default path. Dual A stays **OFF** unless you are running a Dual A session.

## What each file is for

| File                                                                                      | Use in Paper A                                                                           |
|-------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------|
| Research / Share CSV<br>(device_type, iso8601_timestamp, red, ir, ..., spo2_percent, ...) | Wear success, SpO <sub>2</sub> QC, Arm P sketches                                        |
| Clinical Temporalis <b>PDF</b>                                                            | Hypnogram + bands + bout tables (Arm P / overnight)                                      |
| Memory-flush CSVs (app Application Support; hourly)                                       | Rebuild long nights if PDF sample stream thin — John can pull from phone Files if needed |
| Dual A Mac pack / iOS three-file Share                                                    | Methods Dual A + SpO <sub>2</sub> ∇EMG nest ( <b>not</b> AHI)                            |
| CRF                                                                                       | Comfort, setup time, MAD, AcuPebble, adverse events                                      |

## Partner-only (not from Oralable app)

| Item                                     | Owner             |
|------------------------------------------|-------------------|
| CRF complete                             | Pedro / Ed / self |
| MAD type · VDO · advancement             | Pedro             |
| AcuPebble AHI/ODI (nest, do not replace) | Pedro             |
| Bruxoff notes (if used)                  | Pedro             |
| Beacon ethics / consent ID               | Pedro (+ Ed)      |

### Quick checks before you send

---

- Kit ID + date in filename
- FW **1.0.70** / app **4.3.3** on CRF
- Arm named (Core / Arm P /  $\geq 6$  h / Dual A)
- Arm P or overnight → **PDF** attached (not only short continuous CSV)
- Dual A → three files or Mac concordance folder
- No AHI language in filenames or notes for Oralable metrics

### John receives → builds

---

Wear % · usability · SpO<sub>2</sub> QC · Arm P within-subject oxygen sketches · Dual A precursor plots ·  $\geq 6$  h hypnogram when available (replaces FIG-CO-025 layout exemplar).

---

Ops companion to PAPER\_A\_FEASIBILITY\_PROTOCOL §6. Not ethics approval. Source: PAPER\_A\_DATA\_HANDOFF\_SOP.md